


Model DYCP-31BN



Feature


- Made of high quality polycarbonate, exquisite and durable;
- Transparent top lid, easy for observation;
- Auto-switch-off when the lid is opened;
- Removable electrodes, easy to maintain and clean;
- Easy and simple to use;
- Save buffer solution;
- With fluorescent band on the gel tray makes it convenient to add the samples and observe the gel;
- Special comb design, one comb has two different teeth.

Application

Apply to identify, separate, prepare DNA, and measure its molecular weight.

Specification

Dimension (LxWxH)	197×96×100mm
Gel Size (LxW)	60×60mm
Comb	6 wells, 8 wells and 11 wells
Comb Thickness	1.0mm and 1.5mm
Number of Samples	8-11
Buffer Volume	120 ml
Weight	0.5kg



Recommended Power Supply
DYY-6C/6D/8C/10C