


Model DYCP-31CN



Feature

- Made of high quality polycarbonate, exquisite and durable;
- Transparent top lid, easy for observation;
- Special gel casting base for easy and fast gel casting;
- Auto-switch-off when the lid is opened;
- Removable electrodes, easy to maintain and clean;
- Easy and simple to use;
- With black band and fluorescent mark on the gel tray makes it convenient to add the samples and observe the gel;
- Special comb design, one comb has two different teeth.

Application

Apply to identify, separate, prepare DNA, and measure its molecular weight.

Specification

Dimension (LxWxH)	270×115×110mm
Gel Size (LxW)	100×70mm
Comb	8 wells and 15 wells
Comb Thickness	1.0mm and 1.5mm
Number of Samples	8-15
Buffer Volume	260ml
Weight	1.5kg



Recommended Power Supply
DY-6C/6D/8C/10C