


Model DYCP-32B



Feature

- Made of high quality polycarbonate, exquisite and durable;
- Transparent top lid, easy for observation;
- Special gel casting base for easy and fast gel casting;
- Auto-switch-off when the lid is opened;
- Suitable for 12-channel pipette use;
- Removable electrodes, easy to maintain and clean;
- Easy and simple to use;
- High temperature resistant, not easy deformation;
- Save buffer solution;
- With black band and fluorescent mark on the gel tray makes it convenient to add the samples and observe the gel;
- One gel casting base can cast different sizes of gel.

Application

Apply to identify, separate, prepare DNA, and measure its molecular weight.

Specification

Dimension (LxWxH)	365×165×130mm
Gel Size (LxW)	200×130mm, 150×130mm
Comb	14 wells, 26 wells and 18 wells
Comb Thickness	1.0mm and 1.5mm
Number of Samples	4-208
Buffer Volume	800ml
Weight	1.5kg

Recommended Power Supply
DYY-6C/6D/8C