


Model DYCP-38C



Feature

- Delicate appearance;
- The injection molded body, no leakage;
- Three piece of electrodes with platinum wire;
- The level of the tank body can be adjusted;
- Adjusting sticks for different sizes of paper electrophoresis or cellulose acetate membrane electrophoresis experiments.

Application

For paper and slide electrophoresis, cellulose acetate membrane electrophoresis.

Specification

Dimension (LxWxH)	370×270×110mm
Gel Size (LxW)	70 or 90×250mm (Dual-row)
Buffer Volume	1000ml
Weight	2.0kg




Recommended Power Supply
DYY-6C/6D/8C/10C

Cellulose Acetate Membrane

As the necessary product for DYCP-38C electrophoresis cell, we supply cellulose acetate membrane as the following specifications:

Description	Specification	Packing
Soft cellulose acetate membrane (Easy to wet and operate)	70 × 90	50pcs/case
	20 × 80	50pcs/case
	120 × 80	50pcs/case



Note: We can supply customized size of cellulose acetate membrane as your requirements.