

Model DYCP-40E



Feature

- Special design and selected material;
- Semi-dry blotting, no need buffer solution;
- Fast transfer speed and high transfer efficiency;
- Safe plug technique, all the exposed parts are insulated;
- The transfer bands are very clear.

Application

For fast transferring the protein molecule from the gel to membrane like nitrocellulose membrane.

Specification

Dimension (LxWxH)	270×250×104mm
Blotting Area (LxW)	200×200mm
Continuous Working Time	≥24 hours
Weight	5.0kg

Recommended Power Supply
DYY-6C/7C