

Model DYCZ-20A



Feature

- Easy to cast the gel;
- Transparent, no obstructions to visibility;
- Unique design of heat dissipation, keep temperature balance;
- Simple and easy installation of the tank;
- Easy to make gel with the gel filling device;
- Neat and clear electrophoresis bands can be obtained.

Application

For DNA sequencing and DNA fingerprinting analysis, differential display, AFLP or SSCP research.

Specification

Dimension (LxWxH)	240×210×655mm
Gel Size (LxW)	580×170mm
Comb	32 wells (shark teeth) 26 wells (Great Wall teeth)
Comb Thickness	0.4mm
Number of Samples	52-64
Buffer Volume	850ml
Weight	10.5kg

Recommended Power Supply
DYY-10C/12C/12