

Model DYCZ-20C



Feature

- Easy to cast the gel;
- Stable performance;
- Unique design of heat dissipation, keep temperature balance;
- Distinctive outlet for buffer solution;
- Easy to make gel with the gel filling device;
- Neat and clear electrophoresis bands;
- Clear signs on the glass to ensure correct operation.




Application

For DNA sequencing analysis and DNA fingerprinting analysis, differential display and SSCP research.

Specification

Dimension (LxWxH)	410×210×430mm
Gel Size (LxW)	300×340mm
Comb	68 wells (shark teeth)(standard) 100 wells (shark teeth)(optional)
Comb Thickness	0.4mm
Number of Samples	132
Buffer Volume	1500ml
Weight	12kg



Recommended Power Supply
DYY-10C/12C/12