


Model DYCZ-24DN



Feature

- Made of high quality transparent polycarbonate, exquisite and durable, easy for observation;
- With a gel casting in original position, able to cast and run the gel at the same place, simple and convenient to make gels, and save your precious time;
- Special wedge frame design can fix gel room firmly;
- Molded buffer tank equipped pure platinum electrodes;
- Easy and convenient to add samples;
- Able to run one gel or two gels at the same time;
- Save buffer solution;
- Special design of the tank avoid buffer and gel leakage;
- Removable electrodes, easy to maintain and clean;
- Auto-switch-off when the lid is opened.



Application

For protein electrophoresis

Specification

Dimension (LxWxH)	140×100×150mm
Gel Size (LxW)	75×83mm
Comb	10 wells and 15 wells
Comb Thickness	1.0mm and 1.5mm (Standard) 0.75mm (Optional)
Number of Samples	20-30
Buffer Volume	400ml
Weight	1kg



**Recommended Power Supply
DYD-6C/6D/8C/10C/12**