

Model DYCZ-25E



Feature

- Two main body carry 1-4 pieces of gel;
- Glass is optimized design, greatly reduce the breakage possibility;
- Install rubber chamber in the electrophoresis core subject directly, a set of two pieces of glass plate installed respectively;
- Can cast and run 4 gels at the same time, ensure consistent experiment condition;
- Operation requirement is very simple and accurate limit installation design, make high-end product simplification;
- Casting mold with high strength PC material, which is solid durability. Tank is beautiful and transparent, the running status can be shown clearly.

Application

For SDS-PAGE, Protein electrophoresis

Specification

Dimension (LxWxH)		200×175×175mm
Gel Size (LxW)		104×100mm
Comb		13 wells and 21 wells
Comb Thickness		1.0mm and 1.5mm
Number of Samples		52-54
Buffer Volume		850ml (2 gels) 1200ml (4 gels)
Weight		1.5kg

Recommended Power Supply **DYY-6C/6D/8C/10C/12**