

## Model DYCZ-26C



### Feature

- High transparent poly-carbonate plastic injection molded tank;
- With the special gel casting, makes the gel casting easy and reliable;
- Special balance disc keep the gel balance in first dimension;
- Dielectrophoresis can be finished in one day, saving time, lab materials and space;
- No gel and buffer leakage;
- Auto-switch-off when the lid is opened;
- Limit function and accurate operation.

### Application

For 2-DE proteome analysis, need WD-9412A to cooling the second dimension.

### Specification

Dimension (LxWxH)	300 X 160 X 300 mm
Buffer Volume	3500 ml
The 1D Tube Gel Size(DxL)	Φ 2.5 x 90mm or 170mm
The 2D Gel Size(LxW)	200x175mm (dual plates)
Weight	9.5kg



### Recommended Power Supply DYY-10C/12/12C