

## Model DYCZ-40F



### Feature


- Made of high quality transparent polycarbonate with platinum electrodes;
- Fast transfer speed and high transfer efficiency;
- A seamless, injection-molded transparent buffer tank prevents leakage and breakage;
- Customized blue ice pack as cooling unit, which can help the rotor magnetic stirring, better for heat dissipation.
- Compatible with the lid and the buffer tank of DYCZ-25E.

### Application

Used to transfer the protein molecule from the gel to membrane like nitrocellulose membrane in the Western Blot experiment.

### Specification

|                         |               |
|-------------------------|---------------|
| Dimension (LxWxH)       | 200×175×175mm |
| Blotting Area (LxW)     | 104×110mm     |
| Continuous Working Time | ≥24 hours     |
| Buffer Volume           | 1200ml        |
| Weight                  | 1.5kg         |



**Recommended Power Supply**  
**DYY-6C/6D/7C**