



Model DYCP-31DN



Feature

- Made of high quality polycarbonate, exquisite and durable;
- Transparent top lid, easy for observation;
- Special gel casting base for easy and fast gel casting;
- Auto-switch-off when the lid is opened;
- Removable electrodes, easy to maintain and clean;
- Save buffer solution, easy and simple to use;
- With black band and fluorescent mark on the gel tray makes it convenient to add the samples and observe the gel;
- Can run two pieces of gel at the same time;
- One gel casting base can cast different sizes of gel.



**Recommended Power Supply
DY-6C/6D/8C/10C**

Application

Apply to identify, separate, prepare DNA, and measure its molecular weight.

Specification

Dimension (LxWxH)	310×150×120mm
Gel Size (LxW)	60×60mm, 60×120mm, 120×60mm, 120×120mm
Comb	2+3 wells, 6+13 wells, 8+18 wells, and 11+25 wells
Comb Thickness	1.0mm, 1.5mm and 2.0mm
Number of Samples	2-100
Buffer Volume	650ml
Weight	1.0kg